

Name:	Program:	Lab ID :	
Exp No:	Section	Date:	

# **Carey Foster's Bridge (CFB)**

\_\_\_\_\_

Determine the specific resistance of the given wire material.

# MARKS DISTRIBUTION

	Max. Marks	Marks Obtained
Pre - lab preparation	10	
Observation	15	
Calculation	15	
Result	05	
Viva	05	
Total	50	

Staff signature with date





### **CAREY FOSTER'S BRIDGE**

**Aim**: To determine (i) The specific resistance of the material of a wire and (ii) the unknown resistance of the given wire.

**Apparatus Required**: Meter bridge, Leclanché cell, two equal resistances, variable resistance box, unknown resistance, wire, high resistance, switches, galvanometer, jockey.

#### Formula:

- (i) Resistance of the given coil of wire  $X = R + (l_1 l_2) \rho \Omega$  (1)
- (ii) Specific resistance of the given coil of wire  $S = \frac{X\pi r^2}{L} \Omega m$  (2)

Symbol	Description	Unit
X	Unknown resistance	Ω
R	Known value of resistance in the resistance box	Ω
$l_1$ , $l_2$ , $l_3$ and $l_4$	Balancing lengths	m
ρ	Resistivity of the bridge wire	Ω/m
r	Radius of the given coil of wire	m
L	Length of the given wire	m



## **PART I - Determination of ρ (Resistance per unit length)**




Fig 1- R in the left gap and the right gap shorted




Fig 2 – R in the right gap and the left gap shorted

## PART II -Determination of unknown resistance




Fig 3 – Unknown resistance X in the left gap




Fig 4 - Unknown resistance X in the right gap



#### **Principle**:

Let two equal resistances P and Q be connected in gaps  $G_3$  and  $G_4$ . Let a known resistance R be connected in gap  $G_1$  and unknown resistance X in gap  $G_2$ . Let  $I_1$  be the balance length measured from the end A (with known resistance R in the left gap). Therefore balancing length corresponding to X is  $(100 - I_1)$ 

Therefore 
$$\frac{P}{Q} = \frac{R + l_1 \rho + \alpha}{X + (100 - l_1)\rho + \beta}$$
 (3)

Where  $\rho$  = resistance per unit length of the bridge wire and  $\alpha$  and  $\beta$  are the end corrections at the ends A and B. If the experiment is repeated with R and X interchanged in the gaps, and if  $l_2$  is the balancing length measured from the same end A,

$$\frac{P}{Q} = \frac{X + l_2 \rho + \alpha}{R + (100 - l_2)\rho + \beta} \tag{4}$$

From the two equations, we get,

$$X = R + (l_1 - l_2) r_b \quad \Omega$$

Instead of resistance X, let a thick copper strip be connected in the gap  $G_2$ . i.e, X = 0. Let  $l_3$  be the balancing length measured from the end A (with known resistance R in gap  $G_1$ ) With R in gap  $G_2$  and the copper strip connected in  $G_1$ , let  $l_4$  be the balancing length measured from the same end A.

Using the above equation,

$$0 = R + (l_3 \sim l_4)\rho$$
Or, 
$$\rho = \frac{R}{l_3 \sim l_4}$$

Specific resistance,  $S = \frac{\pi r^2 X}{L}$ , where r is the radius and L, the length of the wire.

#### **Procedure:**

### Part I To determine the resistance/unit length $(\rho)$ of the bridge wire.

- a) Connections are made as shown in the Figure 1. P and Q are two equal resistances.
- b) A low range resistance box R is connected in gap G<sub>1</sub> and short circuit to gap G<sub>2</sub> using a thick copper strip of Zero resistance.
  - c) A small resistance (order of  $0.1\Omega$  ) is introduced in R. The balancing length  $l_3$  from the end A is noted.



- d) The resistance R and the copper strip are interchanged in the gaps. (Refer Fig 2). The balancing length  $l_4$  from the end A is again determined.
  - e) Hence  $\rho = \frac{R}{l_3 \sim l_4}$  is calculated. The experiment is repeated for different values of R.
  - f) The mean value of  $\rho$  is calculated.

### Part II To determine the unknown resistance (X):

- a) Instead of copper strip, the coil of unknown resistance of length 1 is connected in gap  $G_1$  (refer Fig 3)
- b) The balancing lengths  $l_1$  and  $l_2$  are determined with R in the right and left extreme gaps respectively.
- c) Knowing  $\rho$ , the value X is calculated using  $X = R + (l_1 l_2)\rho$ .
- d) The experiment is repeated for different values of R and the mean value of X is taken.
- e) Measuring the radius r using screw gauge and using the given the length *l* of the coil, the specific resistance is calculated.

#### **Observations:**

#### (i) Determination of $\rho$

	Balancing Le	ength with R in		n
Resistance, R $(\Omega)$	Left gap $(l_3)$ (m)	Right gap $(l_4)$ (m)	$(l_3 \sim l_4)$ (m)	$\rho = \frac{R}{(l_3 \sim l_4)} \ \Omega/m$

Mean $\rho = \underline{\hspace{1cm}} \Omega / n$
---



### (ii) Determination of X

Resistance, R $(\Omega)$	Balancing Length with R in			$X = R + (l_1 - l_2)\rho$
	Left gap (l <sub>1</sub> ) (m)	Right gap (l <sub>2</sub> ) (m)	$(l_1 - l_2)(\mathbf{m})$	$(\Omega)$

Mean $X =$	$\Omega$
	<del>-</del> -

- (iii) Radius of the wire (r) = \_\_\_\_\_  $x10^{-3}$  m
- (iv) Length of the wire (L) = \_\_\_\_\_m

Specific resistance 
$$S = \frac{\pi r^2 X}{L} = \Omega m$$

#### **Result:**

- (i) The unknown resistance of the given coil of wire (X) =  $\Omega$
- (ii) Specific resistance of the material (S) =  $\Omega$  m

**Signature of Faculty In-charge** 

## Sample viva questions:

- 1. What is the principle of Carey-Foster's bridge?
- 2. What is specific resistance?
- 3. When is Carey-Foster's Bridge more sensitive?



### WORKSHEET